

Soft Lens Solutions

(see note at bottom of page for abbreviations.)

Daily Cleaners

Manufacturer	Name	Preservative
Bausch + Lomb	SENSITIVE EYES Daily Cleaner	sorbic acid .25% and EDTA .5%
Lobob	SOF/PRO2 Extra-Strength Daily Cleaner	preservative-free
Optikem International	SERINE Extra-Strength Daily Cleaner	preservative-free
Orion Vision Group	Lens Fresh Extra Strength	contains isopropyl alcohol 15.7%
www.miraflow.com	MiraFlow Extra-Strength Daily Cleaner	preservative-free

Note: ALCON'S OPTI-FREE Daily cleaner has been discontinued.

Daily Protein Removers

Manufacturer	Name	Preservative
Alcon Vision LLC	OPTI-FREE Replenish Rewetting Drops	Polyquad .001%
Alcon Vision LLC	CLERZ PLUS Lens Drops	lens lubricant also fights protein buildup during lens wear.
Alcon Vision LLC	OPTI-FREE SUPRACLENS DAILY PROTEIN REMOVER	Not for use directly in eye; preservative free; contains pancreatin

Saline Solutions

Non-thimerosal Preserved

Manufacturer	Name	Preservative
Amcon, Inc.	GOOD SENSE Saline Solution	sorbic acid 0.1% & EDTA
Bausch + Lomb	SENSITIVE EYES PLUS Saline Solution	polyaminopropyl biguanide .00003% & EDTA .025%
K-Mart	AMERICAN FARE Saline Solution	polyaminopropyl guanide .00003% & EDTA .025%
Target	Saline Solution for Sensitive Eyes	polyaminopropyl guanide .00003% & EDTA .025%
Walmart	EQUATE Sterile Saline Solution	sorbic acid 0.1% & EDTA .025%

Preservative Free

Manufacturer	Name	Description
Church & Dwight	Simply Saline	12 oz. aerosol; buffered
Menicon	LacriPure	sterile non-preserved saline in 5 ml unit-dose vials
Purilens/Lifestyle Co. Inc.	PURILENS SOLUTION	preservative-free buffered saline in 4 oz. bottles

Lens Lubricants

Manufacturer	Name	Preservative
Advanced Vision Research	THERATEARS Contact Lens Comfort Drops	Sodium perborate
Alcon Vision LLC	OPTI-FREE PUREMOIST Rewetting Drops	Polyquad .001% & edelate disodium 0.05%
Alcon Vision LLC	OPTI-FREE REPLENISH Rewetting Drops	Polyquad .001%
Alcon Vision LLC	OPTI-FREE EXPRESS Rewetting Drops	Polyquad .001%
Allergan	REFRESH RELIEVA for Contacts	preservative-free
Allergan	REFRESH CONTACTS	preservative-free; single use ampules
Allergan	REFRESH CONTACTS	.005% Purite; in .4 fluid oz. bottle
Amcon	VIVA DROPS	sorbic acid 0.25%; EDTA 0.10%
Bausch + Lomb	SENSITIVE EYES Drops	sorbic acid .1%; EDTA
Bausch + Lomb	RENU MULTIPLUS Lubricating & Rewetting Drops	sorbic acid .1% & EDTA .1%
Bausch + Lomb	RENU Rewetting Drops	sorbic acid .1% & EDTA
Eyevance Pharmaceuticals	FRESHCOAT Lubricant Eye Drops	preservative-free 10 ml & single-use vials
Johnson & Johnson Surgical Vision, Inc.	blink Contacts	Ocupure (oxychloro complex) .005%
Johnson & Johnson Surgical Vision, Inc.	Blink-N-Clean Lens Drops	polyhexamethylene biguanide .0001%
K-Mart	AMERICAN FARE Rewetting Drops	sorbic acid .15% & EDTA .1%
Lens Comfort, Inc.	LENS COMFORT Lens Lubricant	N/A
Optics Laboratory	EYE-CEPT Rewetting Drops	preservative-free single-use vials
Target	Lubricating & Rewetting Drops	hexadienoic acid .1% & EDTA .1%
Walmart	EQUATE Sterile Lubricating Drops	sorbic acid .15% & EDTA .1%

Not all available private label products in categories presented on this page are listed due to their prevalence in the market. The ones listed here are presented as samples of their availability.

Other Enzymatic Cleaners

Manufacturer	Name	Description
Johnson & Johnson Surgical Vision, Inc.	ULTRAZYME Enzymatic Cleaner	tablets with subtilisin A for soln. in neutralization of hydrogen peroxide (NOT for thermal or chem. disinfection); safely used only when mixed with sterile saline or other solutions as directed

NOTE: Because the AOA and the AAO strongly recommend that salt tablets NOT be used in cleaning or storing CLs, we do not list salt tablets as contact lens saline products in TYLER'S QUARTERLY.

Oxidation: Solutions/Systems

Manufacturer	Name	Description
Alcon Vision LLC	CLEAR CARE PLUS Cleaning & Disinfection Soln. (No rubbing required for soft lenses)	3% H2O2; 79% NaCl; for use with all lenses; contains catalytic disc for autoneutralization & cup; w/ HydraGlyde Moisture Matrix; red-tipped.
Alcon Vision LLC	CLEAR CARE Cleaning & Disinfection Soln. (No rubbing required for soft lenses)	3% H2O2; 79% NaCl; for use with all lenses; contains catalytic disc for autoneutralization & cup; red-tipped.
CooperVision	REFINE ONE-STEP Cleaning & Disinfection Solution	3% H2O2; one-step preservative-free cleaning & disinfecting solution; w/ lubricating agent for enhanced comfort.
Johnson & Johnson Surgical Vision, Inc.	OXYSEPT ULTRACARE Disinfecting/Neutralizer Soln.	one step disinfection, 3% H2O2; delayed-release catalase (red-tipped); with vitamin B12 color reminder (pink).

UV Subsonic Disinfection

Manufacturer	Name	Description
Purilens/Lifestyle Co. Inc.	PURILENS UV Disinfectant System	cleaning unit uses UV-C radiation to disinfect lenses.

Note: AVT no longer sells Natural lens RDS combination solution.

Note: Alcon has discontinued the Clear Care Rinse & Go saline solution product.

Note: Previously listed AMO Products were purchased by Johnson & Johnson Surgical Vision, Inc.

Combination Solutions

Manufacturer	Name	Preservative
Alcon Vision LLC	OPTI-FREE PURE MOIST MPDS	retains moisture w/ HydraGlyde; (for C, R, D, S, P); contains Polyquad & Aldox
Alcon Vision LLC	OPTI-FREE REPLENISH MPDS (No rubbing required for soft lenses)	retains moisture w/ TearGlyde; (for C, R, D, S, P) contains Polyquad .001% & Aldox.
Alcon Vision LLC	OPTI-FREE EXPRESS MPDS LASTING COMFORT NO RUB FORMULA	Polyquad .001% & Aldox; (P); contains Tetronic 1304 - hydrating agent; "no rubbing needed"
Bausch + Lomb	BioTrue Multi-Purpose Solution	polyaminopropyl biguanide .00013% & polyquaterium .0001% (for C, R, D, P, S)
Bausch + Lomb	renu sensitive	DYMED .00005% for C, R & D
Bausch + Lomb	renu fresh	DYMED .0001%; (for C, R, D, P; & S up to 30 days.)
Clear Conscience, LLC	CLEAR CONSCIENCE Multipurpose Soln.	(for C, R, D, S); polyhexanide .0001%
CooperVision	LITE SOLUTION Multipurpose Soln.	(for C, R, D, S, W, P)
Johnson & Johnson Surgical Vision, Inc.	ACUVUE REVITALENS MPDS	(for C, R, S, D, P)
Johnson & Johnson Surgical Vision, Inc.	COMPLETE Multi-Purpose Solution - Easy Rub	polyhexamethylene biguanide & EDTA (for C, R, S, D, P)
Johnson & Johnson Surgical Vision, Inc.	REVITALENS OCUTEC Disinfecting Soln.	alexadine dihydrochloride .00016% & polyquaterium-1 .0003% (for C, R, S, D, P)
K-Mart	AMERICAN FARE NO Rub Multipurpose Soln.	(for C, R, D, S, P); polyaminopropyl biguanide
K-Mart	AMERICAN FARE NO Rub Multipurpose Soln.	(for C, R, D, S, P); polyaminopropyl biguanide
Tangible Science	Tangible Clean Soln.	(for C & D)
Walmart	EQUATE Sterile Multipurpose Soln.	(for C, R, D, S, P); polyaminopropyl biguanide

NOTE: C - cleaning; R - rinsing; D - disinfecting; S - storing; W - wetting; P - protein removal